NeoGate

[en] Instructions for Use

Flexible aid for full mouth access

English

Models: Regular, Small, X-Small

Intend Use: It is intended to allow a dental professional to retract the lips/cheeks of a patient during an orthodontic or restorative dental procedure for easy access to the field of operation. It is designed to be placed in and around the lips/cheeks. **For single use only.**

Cautions:

1. Please check the retractor for damage or aging before each use.
2. The product is not intended to be reprocessed or reused. The material is not suitable for sterilization. Cleaning and disinfecting agents may have a detrimental effect on the material and lead to undesirable reactions in the patient.

Instruction for use:

1. Evaluate the patient and confirm the type of cheek retractor;
2. Fill the cheek retractor into the mouth to open the mouth;
3. After the treatment , take out of the cheek retractor.

Storage: No exposure and no mixing with corrosive items.
Shelf life: 5 years

Safety information:
In the case of serious incidents related to the product, please contact
SUNGO Cert Gmbh
Address: Harffstr. 47,40591 Düsseldorf, Germany
E-mail:de.rep@sungoglobal.com and your responsible competent authority.

Disposal information:
Remaining stocks must be disposed of according to the corresponding national legal requirements.

Manufacturer Information:

Company:
TANGSHAN ZHENGTONG EXHIBITION CO., LTD.
Address: No.24 Chaoyang Road, Lubei District,
Tangshan City, Hebei Province, China
E-Mail: info@ztdental.com

European Union Authorized Representative

EC REP

Company: SUNGO Cert Gmbh
Address: Harffstr. 47,40591 Düsseldorf, Germany
E-mail:de.rep@sungoglobal.com

Packaging Symbols:

CE	CE symbol
	Manufacturer Name Address
	Manufacture Date
	Use by date
	Distributor
LATEX	Latex Free
+2 +28	Storage temperature
EC REP	Name and Address of European Union Representative
	Recyclable
	See Instruction
MD	Medical Device
2	Do not re-use
	Avoid Sunlight

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Hersteller:

TANGSHAN ZHENGTONG EXHIBITION CO., LTD.
No.24 Chaoyang Road, Lubei District, Tangshan City,
Hebei Province, China
E-Mail: info@ztdental.com

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Produttore:

TANGSHAN ZHENGTONG EXHIBITION CO., LTD.
No.24 Chaoyang Road, Lubei District, Tangshan City,
Hebei Province, Cina
E-Mail: info@ztdental.com

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Proizvajalec:

TANGSHAN ZHENGTONG EXHIBITION CO., LTD.
No.24 Chaoyang Road, Lubei District, Tangshan City,
Hebei Province, Kitajska
E-Mail: info@ztdental.com

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Proizvođač:

TANGSHAN ZHENGTONG EXHIBITION CO., LTD.
No.24 Chaoyang Road, Lubei District, Tangshan City,
Hebei Province, Kina
E-Mail: info@ztdental.com

Ovlašteni predstavnik u Europskoj uniji:

EC REP

SUNGO Cert Gmbh
Harffstr. 47,40591 Düsseldorf, Njemačka
E-mail:de.rep@sungoglobal.com

Simboli na pakiranju:

CE	CE simbol
	Naziv i adresa proizvođača
	Datum proizvodnje
	Rok uporabe / datum isteka valjanosti
	Distributer
LATEX	Bez lateksa
+2 +28	Temperatura skladištenja
EC REP	Naziv i adresa ovlaštenog predstavnika u Europskoj uniji
	Može se reciklirati
	Pogledajte upute za uporabu
MD	Medicinski proizvod
2	Ne koristiti ponovno
	Zaštititi od sunčeve svjetlosti

Rx ONLY

CE

Date information prepared:
2026-07-15 / Rev. 1